UNDERWAX
UNDERWAX
UNDERWAX
UNDERWAX
UNDERWAX
UNDERWAX
UNDERWAX
UNDERWAX
UNDERWAX
UNDERWAX

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Master model to carry out the rigid coping: the best solution to obtain a structure of light curing material that completely burns out and can be milled

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Rigid coping carried out with Underwax Nobil Metal

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Underwax is the innovative pattern resin that, thanks to its high stability, prevents from the serious drawbacks of distortion noticed following the conventional technique (wax-up).

Underwax belongs to STRATEGIC-LINE, the protocol that considers all the process technical phases, from the development of the model to the metal-ceramic interface, and which, used as a reference, guarantees a continuity of the result.

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The modelled piece carried out with UNDERWAX must be considerably reduced to have less resin inside the investment mass. It is recommended not to exceed a thickness of 1 mm; should a larger thickness be necessary, proceed with the polymerisation in layers (1 mm at a time). If necessary, interrupt the resin copings by carrying out cuttings (vestibular lingual), which will then be united with wax: in this way, a correct expansion is obtained during setting of the investment and of the entire modelled piece. We recommend using a traditional investment following the manufacturer's suggested firing cycles to allow a complete calcination of the resin components.

«Speed» investment: speed investments, like Proper Vest Fast NF by Nobil Metal, are recommended for minimal masses and strictly observing the manufacturer's directions. Using a speed investment, carry out a slow, conventional preheating cycle.

Caution: photopolymerisable systems with UV light emission do not polymerise the product. Upon completion of the operation, carefully re-close the syringe and store away from light and heat sources.

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Element waxed-up on a rigid coping of Underwax Nobil Metal

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
Frameworks carried out with Underwax and converted into metal

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Thickness check and work ready for spruing

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Pattern resin; wax-up viscous liquid containing activator

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• SOSTANZE PERICOLOSE – DANGEROUS COMPONENTS

Nessuna – None

• PACKAGING E CONTENUTO – PACKAGING AND CONTENT

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We would like to thank Mr ANGELO GALLOCCHIO for the iconographic material.

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See the material safety data sheet, please follow the suggestions on storage and check the expiry date.

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NOBIL-METAL® S.p.A.