

INCISIVI SUPERIORI
 UPPER INCISOR
 INCISIVES SUPERIEURES
 OBERE SCHNEIDEZÄHNE
 INCISIVOS SUPERIORES



ACRY STAR

E 1		
(1) 43,3 (2) 8,7 (3) 10,4 (4) 8,1		

E 2		
(1) 42,7 (2) 9,2 (3) 11,0 (4) 8,1		

E 3		
(1) 45,0 (2) 10,1 (3) 12,1 (4) 8,4		

E 4		
(1) 48,8 (2) 10,5 (3) 12,7 (4) 8,8		

E 5		
(1) 52,0 (2) 13,0 (3) 14,6 (4) 9,2		

G 1		
(1) 40,8 (2) 8,0 (3) 9,7 (4) 7,3		

G 2		
(1) 44,5 (2) 9,2 (3) 10,5 (4) 8,2		

G 3		
(1) 44,8 (2) 9,8 (3) 11,6 (4) 8,2		

G 4		
(1) 48,3 (2) 9,3 (3) 11,8 (4) 8,7		

F 1		
(1) 38,8 (2) 8,8 (3) 11,1 (4) 7,2		

F 2		
(1) 42,1 (2) 8,9 (3) 11,6 (4) 7,7		

F 3		
(1) 42,2 (2) 9,8 (3) 12,4 (4) 7,9		


F 4		
(1) 48,3 (2) 10,4 (3) 12,5 (4) 8,6		

H 1		
(1) 40,6 (2) 8,9 (3) 10,4 (4) 7,7		

H 2		
(1) 42,0 (2) 9,1 (3) 11,4 (4) 7,5		

H 3		
(1) 42,6 (2) 8,8 (3) 10,4 (4) 7,7		

H 4		
(1) 43,8 (2) 9,0 (3) 11,4 (4) 8,0		




INCISIVI INFERIORI
 LOWER INCISORS
 INCISIVES INFÉRIEURES
 UNTERE SCHNEIDEZÄHNE
 INCISIVOS INFERIORES

COMBINATED ARTICULATIONS		
ANTERIOR		POSTERIOR
Upper	Lower	
E1	K5 K3	M2
E2	K6 K5	M4
E3	K6	M4
E4	K9	M5
E5	K10	M5
F1	K1	M1
F2	K2	M2
F3	K5 K2	M3
F4	K10 K9	M5
G1	K2 K1	M1
G2	K6 K7	M2
G3	K6 K4	M4
G4	K8 K1	M5
H1	K7 K1	M1
H2	K5 K2	M2
H3	K3 K1	M3
H4	K4 K3 K5	M4

AVAILABLE IN THE FOLLOWING COLORS	
11 13 WHITE 15	31 BROWN 32 BROWN
20 GREY 21 GREY	39 GREYISH 40 GREYISH 41 BROWN 42 BROWN
17 YELLOW 22 23 25 PINK 27	

K 1 (1) 32,3 (2) 8,1 (3) 9,7 (4) 4,4		
K 2 (1) 32,7 (2) 10,0 (3) 11,4 (4) 4,5		
K 3 (1) 33,8 (2) 8,5 (3) 9,7 (4) 4,6		
K 4 (1) 34,6 (2) 9,8 (3) 11,3 (4) 5,0		
K 5 (1) 34,0 (2) 9,0 (3) 10,9 (4) 5,1		
K 6 (1) 35,0 (2) 9,8 (3) 11,0 (4) 5,3		
K 7 (1) 34,8 (2) 8,4 (3) 9,6 (4) 5,0		
K 8 (1) 40,0 (2) 8,2 (3) 10,1 (4) 6,1		
K 9 (1) 38,5 (2) 9,5 (3) 11,2 (4) 6,5		
K 10 (1) 40,0 (2) 12,2 (3) 14,0 (4) 5,7		

DIATORICI SUPERIORI E INFERIORI
 UPPER AND LOWER POSTERIOR
 MOLAIRES SUPERIEURES ET INFERIEURES
 OBERE UND UNTERE BACKENZÄHNE
 MOLARES SUPERIORES Y INFERIORES



ACRY STAR

ACRY STAR teeth comply with ISO 3336


M 1 (1) 59,2 - (2) 61,7



M 2 (1) 59,4 - (2) 62,8




M 3 (1) 62,3 - (2) 67,2



M 4 (1) 62,6 - (2) 67,3



M 5 (1) 67,0 - (2) 71,5




RUTHINIUM


Dental manufacturing

Via Cà Mignola Nuova, 1699
 Badia Polesine (Rovigo) Italy
 Tel. (0425) 51628-590587
 Telex 430841 RUTBAI I, - Fax (0425) 590156