



## CZR FULL KIT

CZR (CERABIEN ZR) is specially developed porcelain to make all ceramic crowns and bridges, using Zirconia frameworks. It does not cause any allergies and has an advantage not to show black lines at all in the cervical area. The thermal expansion of CZR is totally different from that of other porcelain.

Therefore, mixing or using with other porcelain is not recommended.

Coefficient of thermal expansion  
(50—500°C 10<sup>-6</sup> K<sup>-1</sup>)

EX-3	Cerabien	CZR
12.4	6.8	9.1

SHADE BASE	10g each	18 colors
BODY		18 colors
ENAMEL		3 colors
TRANSLUCENT		4 colors
CERVICAL		4 colors
MODIFIER		14 colors
ADD-ON		2 colors
COLOR COMBINATION TABLE		1 piece
TECHNICAL INSTRUCTIONS		1 copy



## CZR NP SHADE KIT

CZR NP SHADE KIT is developed to create two Noritake original shades (NP1.5 & NP2.5). These shades are often seen in natural teeth.

NP-1.5 : The chroma is between A1 and A2 shades - the hue corresponds to the slightly reddish shades.

NP-2.5 : The chroma is between A2 and A3 shades - the hue corresponds to the slightly reddish shades.

SHADE BASE	10g each	2 colors
OPACIOUS BODY		2 colors
BODY		2 colors
MARGIN		2 colors



## CZR MARGIN KIT

MARGIN PORCELAIN is developed exclusively for the combination application with Alumina framework. In addition to ordinary usage, it can be used to correct a chip in the Alumina framework. The use of MAGIC FORMER as forming liquid is highly recommended.

BASIC COLOR	10g each	14 colors
SPECIAL COLOR		3 colors
DILUTION POWDER (MDL)		1 color
RETOUCHING POWDER (MRP)		1 color



## CZR OPACIOUS BODY KIT

OPACIOUS BODY can be used to achieve delicate colors and to control translucency.

BASIC COLOR	10g each	16 colors
SPECIAL COLOR		3 colors



## CZR LUSTER KIT

LUSTER PORCELAIN reproduces natural enamel in both surface features and color.

- LUSTER PORCELAIN reproduces the fine surface structure and luster of natural teeth.
- A unique combination of fine surface particles produces a selective reflection of light that results in the same opalescence seen in natural teeth.
- LUSTER PORCELAIN has transparent, bright and vivid colors, therefore darkening at the incisal edge or at the occlusal surface will not occur.
- Noritake has thoroughly studied color changes in natural teeth caused by aging. LUSTER PORCELAIN features a complete line of colors consistent with these changes.

BASIC COLOR	10g each	4 colors
SPECIAL COLOR		4 colors



## CZR INTERNAL STAIN KIT

Using INTERNAL STAIN, it is possible to replicate many of the most intricate three-dimensional patterns of shade variation found in natural teeth.

1. WHITE	7. CERVICAL 2	13. A+	3g each
2. INCISAL BLUE 1	8. CERVICAL 3	14. B+	
3. INCISAL BLUE 2	9. EARTH BROWN	15. C+	
4. MAVEION ORANGE 1	10. REDDISH BROWN	16. D+	
5. MAVEION ORANGE 2	11. SALMON PINK		
6. CERVICAL 1	12. RED		
BRIGHT (DILUTION)		3g	1 color
IS LIQUID		10ml	1 bottle
IS COLOR GUIDE			1 piece
TECHNICAL INSTRUCTIONS			1 copy



## CZR EXTERNAL STAIN KIT

EXTERNAL STAIN reproduces the great variety of color shades observed on the surface of natural teeth.

1. PURE WHITE	8. ORANGE 1	15. SALMON PINK	3g each
2. GRAY	9. ORANGE 2	16. PINK	
3. BLACK	10. CERVICAL 1	17. RED	
4. BLUE	11. CERVICAL 2	18. A+	
5. GREEN 1	12. CERVICAL 3	19. B+	
6. GREEN 2	13. EARTH BROWN	20. C+	
7. YELLOW	14. REDDISH BROWN	21. D+	
GLAZE		10g	1 color
ES LIQUID		10ml	1 bottle
ES COLOR GUIDE			1 piece
TECHNICAL INSTRUCTIONS			1 copy

## CZR FORMING LIQUID

Specially formulated for CZR.

100ml

## IS LIQUID

Specially formulated for Internal Stain. It can be used for CZR PRESS Shade Base Stain & CZR PRESS Glaze powder.

10ml

## MEISTER LIQUID

This liquid is slow-drying forming liquid. Keep porcelain moist, which makes applying Shade Base or other powders easily.

100ml

## ES LIQUID

Keeps EXTERNAL STAIN moist, which makes staining work easy.

10ml