



**CZR PRESS SHADE BASE STAIN KIT**

CZR PRESS SHADE BASE STAIN KIT is a kind of the stain to make the basic shade color of zirconia framework before pressing CZR PRESS INGOT. The basic shade color can be recreated thinly, transparently and naturally because of the stain.

BASIC COLOR	3g	20 colors
IS LIQUID	10ml	1 bottle
SB STAIN COLOR GUIDE		1 piece
TECHNICAL INSTRUCTIONS		1 copy



**CZR PRESS INGOT-L KIT / INGOT-H KIT**

CZR PRESS INGOT-L KIT is a low translucent type of press ceramic for zirconia framework. It is developed to recreate the margin area precisely by pressing INGOT-L.

INGOT-L is a press ceramic using with the layering technique. After INGOT-L is pressed, the crown can be made by using CZR TRANSLUCENT PORCELAIN or CZR LUSTER PORCELAIN.

Can make a single anterior crown without frameworks using with CZR PRESS LF porcelain.

BASIC COLOR	10 g	20 colors
PLUNGER		1 set
TECHNICAL INSTRUCTIONS		1 copy

5g Ingot is available as individual item.



CZR PRESS INGOT-H KIT is a high translucent type of press ceramic for zirconia framework. It is developed to recreate the margin area precisely by pressing INGOT-H.

INGOT-H is a press ceramic using with the stain technique. After INGOT-H is pressed, the crown can be made by using CZR EXTERNAL STAIN or CZR PRESS GLAZE POWDER.

CZR PRESS LF porcelain can be used after pressing also.

BASIC COLOR	10 g	20 colors
PLUNGER		1 set
TECHNICAL INSTRUCTIONS		1 copy

5g Ingot is available as individual item.



**CZR PRESS TRIAL-L KIT / TRIAL-H KIT**

CZR PRESS TRIAL-L KIT is a necessary kit for using INGOT-L, including SHADE BASE STAIN, INGOT-L, CZR PORCELAIN, and PLUNGER. It contains 5 shade colors, which are A1, A2, A3, NW0, and NW0.5. The beginning user of CZR PRESS should also prepare for NORITAKE PRESS INVESTMENT, RING / FORMER, IS LIQUID, and CZR FORMING LIQUID.

INGOT-L	10g	5 colors
SB STAIN	3g	5 colors
CZR PORCELAIN POWDER	10g	8 colors
PLUNGER		1 set
SB STAIN COLOR GUIDE		1 piece
CZR PORCELAIN INSTRUCTION		1 copy
TECHNICAL INSTRUCTIONS		1 copy



CZR PRESS TRIAL-H KIT is a necessary kit for using INGOT-H, including SHADE BASE STAIN, INGOT-H, CZR PRESS GLAZE, CZR ES (A+), and PLUNGER. It contains 5 shade colors, which are A1, A2, A3, NW0, and NW0.5. The beginning user of CZR PRESS should also prepare for NORITAKE PRESS INVESTMENT, RING / FORMER, IS LIQUID, and ES LIQUID.

INGOT-H	10g	5 colors
SB STAIN	3g	5 colors
CZR ES A+	3g	1 color
CZR PRESS GLAZE	10g	1 color
PLUNGER		1 set
SB STAIN COLOR GUIDE		1 piece
CZR PORCELAIN INSTRUCTION		1 copy
TECHNICAL INSTRUCTIONS		1 copy



CZR PRESS INTRO-L KIT is the optimum kit for those who first start to use the CZR PRESS SYSTEM. It contains 5 shade colors including A1, A2, A3, NWO, and NWO.5.



CZR PRESS INTRO-H KIT is the optimum kit for those who first start to use the CZR PRESS SYSTEM. It contains 5 shade colors including A1, A2, A3, NWO, and NWO.5.

**CZR PRESS INTRO-L KIT / INTRO-H KIT**

**INTRO-L KIT**

CZR PRESS TRIAL-L KIT		1 set
RING/FORMER SET 100/200		1 set
PRESS INVESTMENT SET	3 kg	1 set
CZR FORMING LIQUID	10ml	1 bottle
IS LIQUID	10ml	1 bottle

**INTRO-H KIT**

CZR PRESS TRIAL-H KIT		1 set
RING/FORMER SET 100/200		1 set
PRESS INVESTMENT SET	3 kg	1 set
CZR ES LIQUID	10ml	1 bottle
CZR IS LIQUID	10ml	1 bottle



**CZR PRESS LAYERING KIT**

CZR LAYERING KIT is CZR porcelain selected for using CZR PRESS INGOT-L.

BODY	10g each	20 colors
ENAMEL		3 colors
TRANSLUENT		3 colors
MODIFIER		14 colors
ADD-ON		2 colors
LUSTER		8 colors
MRP		1 color
COLOR TABLE/WITHOUT SB		1 piece
TECHNICAL INSTRUCTIONS		1 copy



**CZR PRESS LF ASSORTMENT KIT**

Noritake CZR PRESS LF is low fusing porcelain to be built up to make an enamel layer after pressing CZR PRESS ingot. By using this porcelain with pressed ceramic without a Zirconia framework, you can make a single crown, a porcelain laminate veneer, an inlay and an onlay. Please use Noritake LF Liquid for this porcelain.

ENAMEL	10g each	3 colors
TRANSLUENT/LUSTER		12 colors
AD-T		1 color
TECHNICAL INSTRUCTIONS		1 copy

**LF LIQUID**

Special Liquid for LF porcelain.

100ml
-------

AD-T will be included.

**CZR PRESS GLAZE POWDER**



This CZR PRESS GLAZE POWDER is necessary to produce glossiness of the crown in the H-ingots staining technique. Be sure to use this glaze powder and refer to our CZR PRESS technical instruction for the glazing temperature. This glaze powder is available individually. **The use of IS Liquid is highly recommended.**

**RING/FORMER SET**



NORITAKE RING/FORMER SET is a RING/FORMER used when pressing CZR PRESS INGOT. 3 sizes are available, 100g, 200g and 300g.

**RING/FORMER SET 100**

RING 100	FORMER 100
----------	------------

**RING/FORMER SET 200**

RING 200	FORMER 200
----------	------------

**RING/FORMER SET 300**

RING 300	FORMER 300
----------	------------

**RING/FORMER SET 100/200**

RING 100	FORMER 100
RING 200	FORMER 200

**NORITAKE PLUNGER**



NORITAKE PLUNGER is a plunger used when pressing CZR PRESS INGOT.

3 pieces
----------

**NORITAKE DISPO PLUNGER**



Noritake DISPO PLUNGER is disposable plunger set. 2 kinds of size are available for each 2g and 5g CZR PRESS INGOT.

φ12mm/H37mm	50 pcs.
φ16.5mm/H42mm	50 pcs.