



NORITAKE SUPER PORCELAIN EX-3 PRESS PASTE OPAQUE BASIC KIT

EX-3 PRESS PASTE OPAQUE BASIC KIT is specially designed for EX-3 PRESS INGOT. It can make the basic shade color of metal framework with a thinner coat before pressing ingot.

EX3-P POnA1	EX3-P POnB1	EX3-P POnC2	EX3-P POnD4
EX3-P POnA2	EX3-P POnB2	EX3-P POnC3	EX3-P PONW0
EX3-P POnA3	EX3-P POnB3	EX3-P POnC4	EX3-P PONW0.5
EX3-P POnA3.5	EX3-P POnB4	EX3-P POnD2	EX3-P PONP1.5
EX3-P POnA4	EX3-P POnC1	EX3-P POnD3	EX3-P PONP2.5
TECHNICAL INSTRUCTIONS			1 copy
PASTE OPAQUE MODIFIER BRUSH			1 piece



NORITAKE SUPER PORCELAIN EX-3 PRESS PASTE OPAQUE MODIFIER KIT

NORITAKE EX-3 PRESS PASTE OPAQUE MODIFIER KIT contains 8 color modifiers which may be incrementally incorporated into EX-3 PRESS Basic Paste Opaques to create custom shades. They may also be applied individually for minor color adjustments & subtle shade variances.

1. WHITE	5. REDDISH BROWN	3g each
2. GRAY	6. PINK	
3. ORANGE	7. BLUE	
4. EARTH BROWN	8. YELLOW	
TECHNICAL INSTRUCTIONS		1 copy
PASTE OPAQUE MODIFIER BRUSH		1 piece



NORITAKE SUPER PORCELAIN EX-3 PRESS INTRO-L KIT / INTRO-H KIT

NORITAKE EX-3 PRESS INTRO-L KIT is an ideal kit for the new user who wishes to try the Ingot-L EX-3 PRESS & a sampling of the EX-3 LF Layering assortment. It includes everything in the EX-3 PRESS INGOT-L Trial kit, plus: Noritake Press Investment, Rings / formers & LF Liquid. This kit is for the technician who desires a 'turn-key' experience.

LnA1	EX3-P POnA1	EX3-P LF E1	EX3-P LF Tx
LnA2	EX3-P POnA2	EX3-P LF E2	EX3-P LF AD-T
LnA3	EX3-P POnA3	EX3-P LF E3	EX3-P LF MRP
LNW0	EX3-P PONW0	EX3-P LF LT1	
LNW0.5	EX3-P PONW0.5	EX3-P LF CE	
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces
RING / FORMER SET			1 set
PRESS INVESTMENT SET			1 set
LF FORMING LIQUID			10ml
PASTE OPAQUE MODIFIER BRUSH			1 piece



NORITAKE EX-3 PRESS INTRO-H KIT is an ideal kit for the new user who wishes to try the Ingot-H EX-3 PRESS & the EX-3 LF External Stain & Glaze. It includes everything in the EX-3 PRESS INGOT-H Trial kit, plus: Noritake Press Investment, Rings / formers & ES Liquid. This kit is for the technician who desires a 'turn-key' experience.

HnA1	HNW0	EX3-P POnA2	EX3-P PONW0.5
HnA2	HNW0.5	EX3-P POnA3	EX3-P ES A+
HnA3	EX3-P POnA1	EX3-P PONW0	EX3-P GLAZE
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces
RING / FORMER SET			1 set
PRESS INVESTMENT SET			1 set
ES LIQUID			10ml
IS LIQUID			10ml
PASTE OPAQUE MODIFIER BRUSH			1 piece



NORITAKE SUPER PORCELAIN EX-3 PRESS INGOT-L KIT / INGOT-H KIT

NORITAKE EX-3 PRESS INGOT-L Kit contains low translucency pressable ceramic ingots. It was developed to recreate both the margin area, as well as the reduced dentinal form of the tooth. The L-ingot is indicated for use in the layering method. After pressing the L-ingot, the enamel layer of the tooth can be created by using EX-3 PRESS LF Translucent & Luster porcelains.

LnA1	LnB1	LnC2	LnD4
LnA2	LnB2	LnC3	LNW0
LnA3	LnB3	LnC4	LNW0.5
LnA3.5	LnB4	LnD2	LNPN1.5
LnA4	LnC1	LnD3	LNPN2.5
TECHNICAL INSTRUCTIONS			1 copy
EX-3 DISPO PLUNGER			3 pieces



NORITAKE EX-3 PRESS INGOT-H Kit contains high translucency pressable ceramic ingots. It was developed to recreate both the margin area, as well as the full-contour form of the tooth. The H ingot is indicated for use in the staining method. After pressing the H ingot, the restoration may be characterized using EX-3 PRESS LF External Stain & EX-3 PRESS Glaze Powder. The EX-3 PRESS LF Layering Assortment may also be used after pressing.

HnA1	HnB1	HnC2	HnD4
HnA2	HnB2	HnC3	HNW0
HnA3	HnB3	HnC4	HNW0.5
HnA3.5	HnB4	HnD2	HNP1.5
HnA4	HnC1	HnD3	HNP2.5
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces



NORITAKE SUPER PORCELAIN EX-3 PRESS TRIAL-L KIT / TRIAL-H KIT

NORITAKE EX-3 PRESS TRIAL-L KIT is an introductory kit which includes EX-3 PRESS Opaque, L-ingots, EX-3 PRESS LF Enamel, Luster & Translucent Porcelains. It contains the following 5 shades: NW0, NW0.5, A1, A2 & A3. The beginner EX-3 PRESS user should also use Noritake PRESS Investment, Ring / former, Paste Opaque Liquid & LF Liquid.

LnA1	EX3-P PONA1	EX3-P LF E1	EX3-P LF Tx
LnA2	EX3-P PONA2	EX3-P LF E2	EX3-P LF ADT
LnA3	EX3-P PONA3	EX3-P LF E3	EX3-P LF MRP
LNW0	EX3-P PONW0	EX3-P LF LT1	
LNW0.5	EX3-P PONW0.5	EX3-P LF CE	
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces



EX-3 PRESS TRIAL-H KIT is an introductory kit which includes EX-3 PRESS Opaque, H ingots, EX-3 PRESS LF External Stain A+ & EX-3 PRESS LF Glaze Powder. It contains the following 5 shades: NW0, NW0.5, A1, A2 & A3. The beginner EX-3 PRESS user should also use Noritake PRESS Investment, Ring / former, Paste Opaque Liquid & ES Liquid.

HnA1	HNW0	EX3-P PONA2	EX3-P PONW0.5
HnA2	HNW0.5	EX3-P PONA3	EX-3 PRESS A+
HnA3	EX3-P PONA1	EX3-P PONW0	EX-3 PRESS GLAZE
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces



NORITAKE SUPER PORCELAIN EX-3 PRESS ESTHETIC WHITE KIT

NORITAKE EX-3 PRESS ESTHETIC WHITE KIT is a supplemental esthetic kit which contains the necessary materials for reproducing the complex bleach shades today's estheticdriven patients present. In addition to brighter & whiter press opaques & pressable ingots, this kit contains two 'must-have,' high-value enamels, which reproduce the enamel appearance many patients exhibit as a result of having bleached their teeth.

INGOTS 2g x 5	EW00, EYW0, EYW0.5, EYW1	
PASTE OPAQUE 10g	POEW, POEYW	
EX3-P LF 10g	SILKY E1, SILKY E2	
TECHNICAL INSTRUCTIONS		1 copy
EX-3 PRESS DISPO PLUNGER		3 pieces



NORITAKE SUPER PORCELAIN EX-3 PRESS LF BODY KIT

EX-3 PRESS LF BODY KIT is an assortment of low fusing body porcelain to be built up after pressing EX-3 PRESS INGOT.

EX3-P LF nA1B	EX3-P LF nB1B	EX3-P LF nC2B	EX3-P LF nD4B
EX3-P LF nA2B	EX3-P LF nB2B	EX3-P LF nC3B	EX3-P LF NP1.sB
EX3-P LF nA3B	EX3-P LF nB3B	EX3-P LF nC4B	EX3-P LF NP2.sB
EX3-P LF nA3.sB	EX3-P LF nB4B	EX3-P LF nD2B	EX3-P LF NWoB
EX3-P LF nA4B	EX3-P LF nC1B	EX3-P LF nD3B	EX3-P LF NWo.sB
TECHNICAL INSTRUCTIONS			1 copy
LF BODY COLOR GUIDE			1 piece



NORITAKE SUPER PORCELAIN EX-3 PRESS LF ASSORTMENT KIT

NORITAKE EX-3 PRESS LF consists of an assortment of low-fusing basic enamels, translucent enamels & opalescent translucent enamels specifically created for layering EX-3 PRESS. These premium low-fusing ceramics enable the ceramist to recreate the intricate enamel canvas, including delicate transparent areas, tinted incisal subtleties and soft opalescent nuances.

EX3-P LF E1	EX3-P LF T1	EX3-P LF LT NATURAL	EX3-P LF ADT
EX3-P LF E2	EX3-P LF T2	EX3-P LF CREAMY ENAMEL	EX3-P LF ADB
EX3-P LF E3	EX3-P LF LT0	EX3-P LF SUN BRIGHT	EX3-P LF MRP
EX3-P LF Tx	EX3-P LF LT1	EX3-P LF INCISAL AUREOLA	
EX3-P LF To	EX3-P LF T-BLUE	EX3-P LF CREAMY WHITE	
TECHNICAL INSTRUCTIONS			1 copy
LF E & T COLOR GUIDE			1 piece



NORITAKE SUPER PORCELAIN EX-3 PRESS LF INTERNAL STAIN KIT

NORITAKE EX-3 PRESS LF INTERNAL STAIN is specifically & exclusively formulated for internal staining of EX-3 PRESS restorations. EX-3 LF Internal Stain facilitate the reproduction of subtle color nuances, high & low value areas & delicate characterizations visible in natural teeth.

WHITE	CERVICAL 2	A+	3g each
INCISAL BLUE 1	CERVICAL 3	B+	
INCISAL BLUE 2	EARTH BROWN	C+	
MAMELON ORANGE 1	REDDISH BROWN	D+	
MAMELON ORANGE 2	SALMON PINK		
CERVICAL 1	RED		
BRIGHT		3g	1 bottle
IS LIQUID		10ml	1 bottle
IS COLOR GUIDE			1 piece
TECHNICAL INSTRUCTIONS			1 copy



NORITAKE SUPER PORCELAIN EX-3 PRESS LF EXTERNAL STAIN KIT

NORITAKE EX-3 PRESS LF EXTERNAL STAIN KIT is specifically & exclusively formulated for external staining of EX-3 PRESS & EX-3 PRESS stand-alone, pressable restorations. EX-3 LF External Stains enable the ceramist to reproduce the complex surface colors, stains & subtle characterizations visible on the surface of natural teeth. Due to their precise thermal-coefficient match & low-fusing firing temp., EX-3 PRESS LF External Stain & EX-3 PRESS LF Glaze are the perfect compliments to EX-3 PRESS & EX-3 PRESS stand-alone, press-able restorations.

PURE WHITE	ORANGE 1	SALMON PINK	3g each
GRAY	ORANGE 2	PINK	
BLACK	CERVICAL 1	RED	
BLUE	CERVICAL 2	A+	
GREEN 1	CERVICAL 3	B+	
GREEN 2	EARTH BROWN	C+	
YELLOW	REDDISH BROWN	D+	
GLAZE		10g	1 color
ES LIQUID		10ml	1 bottle
ES COLOR GUIDE			1 piece
TECHNICAL INSTRUCTIONS			1 copy

NORITAKE SUPER PORCELAIN EX-3 PRESS DISPO PLUNGER

NORITAKE EX-3 PRESS DISPO PLUNGER is disposable plunger set for EX-3 PRESS. 2 kinds of size are available for each 2g & 5g EX-3 PRESS INGOT.