

**IMPIANTI DENTALI**  
**DENTAL IMPLANTS**



## Compatibile Biomet 3i™

Compatible Biomet 3i™

Connessione esagono esterno  
External hexagon connection




	Ø 3	Ø 3,4	Ø 3,8	Ø 4,5	Ø 5,5
H					
8	•	EE340080B	EE380080B	EE450080B	EE550080B
9,5	•	EE340095B	EE380095B	EE450095B	EE550095B
11	•	EE340110B	EE380110B	EE450110B	EE550110B
13	•	EE340130B	EE380130B	EE450130B	EE550130B
15	•	EE340150B	EE380150B	EE450150B	•
Connessione Connection	•	• Ø 3,25	• Ø 3,75-4	• Ø 3,75-4	• Ø 5
	•	• Biomet	• Biomet	• Biomet	• Biomet

## Compatibile Certain®

Compatible Certain®

Connessione esagono interno  
Internal hexagon connection




	Ø 3	Ø 3,4	Ø 3,8	Ø 4,5	Ø 5,5
H					
8	•	EI340080C	EI380080C	EI450080C	EI550080C
9,5	•	EI340095C	EI380095C	EI450095C	EI550095C
11	•	EIE340110C	EI380110C	EI450110C	EI550110C
13	•	EI340130C	EI380130C	EI450130C	EI550130C
15	•	EI340150C	EI380150C	EI450150C	•
Connessione Connection	•	• Ø 3,25	• Ø 3,8	• Ø 4	• Ø 5
	•	• Certain	• Certain	• Certain	• Certain

## Compatibile Xive®

Compatible Xive®


Connessione esagono interno  
Internal hexagon connection



	Ø 3	Ø 3,4	Ø 3,8	Ø 4,5	Ø 5,5
H					
8	EI300080X	EI340080X	EI380080X	EI450080X	EI550080X
9,5	EI300095X	EI340095X	EI380095X	EI450095X	EI550095X
11	EI300110X	EIE340110X	EI380110X	EI450110X	EI550110X
13	EI300130X	EI340130X	EI380130X	EI450130X	EI550130X
15	EI300150X	EI340150X	EI380150X	EI450150X	•
Connessione Connection	• Ø 3	• Ø 3,4	• Ø 3,8	• Ø 4,5	• Ø 5,5
	• Xive	• Xive	• Xive	• Xive	• Xive

### Compatibile Sweden® Compatible Sweden®

Connessione esagono interno  
Internal hexagon connection




		Ø 2,7	Ø 3,4	Ø 3,8	Ø 4,5	Ø 5,5
H	8	•	EI340080S	EI380080S	EI450080S	EI550080S
	9,5	•	EI340095S	EI380095S	EI450095S	EI550095S
	11	•	EI340110S	EI380110S	EI450110S	EI550110S
	13	•	EI340130S	EI380130S	EI450130S	EI550130S
	15	•	EI340150S	EI380150S	EI450150S	•
	Connessione Connection	•		• Ø 3,3	• Ø 3,8	• Ø 4,25
•			■ Sweden	■ Sweden	■ Sweden	■ Sweden

		Ø 2,7	Ø 3,4	Ø 3,8	Ø 4,5	Ø 5,5
H	8	•	EI340080Z	EI380080Z	EI450080Z	EI550080Z
	9,5	•	EI340095Z	EI380095Z	EI450095Z	EI550095Z
	11	•	EIE340110Z	EI380110Z	EI450110Z	EI550110Z
	13	•	EI340130Z	EI380130Z	EI450130Z	EI550130Z
	15	•	EI340150Z	EI380150Z	EI450150Z	•
	Connessione Connection	•		• Ø 3,75	• Ø 3,75	• Ø 4,2
•			■ Zimmer	■ Zimmer	■ Zimmer	■ Zimmer

### Compatibile Zimmer® Compatible Zimmer®

Connessione esagono interno  
Internal hexagon connection



### Standard

Connessione esagono esterno  
External hexagon connection



	-	Ø 3,25	Ø 3,75	Ø 4	Ø 5
<b>H</b>					
5	•	•		EE400050J	EE500050J
6	•	•	EE375060J	EE400060J	EE500060J
7	•	•	EE375070J	EE400070J	EE500070J
8,5	•	EE 375 085	EE375085J	EE400085J	EE500085J
10	•	EE 375 100	EE375100J	EE400100J	EE500100J
11,5	•	EE 375 115	EE375115J	EE400115J	EE500115J
13	•	EE 375 130	EE375130J	EE400130J	EE500130J
15	•	EE 375 150	EE375150J	EE400150J	•
18	•	•	EE375180J	EE400180J	•
<b>Piattaforma Platform</b>	•	• Ø3,4	• Ø4	• Ø4	• Ø5

### Micro Mini


Connessione esagono esterno  
External hexagon connection



	Ø 2,6	Ø 2,9	Ø 3,25	Ø 4	Ø 5
<b>H</b>					
6,5	•	•	•	•	•
7,5	•	•	EE325070D	•	•
8,5	•	•	EE325085D	•	•
10	•	•	EE325100D	•	•
11,5	•	•	EE325115D	•	•
13	•	•	EE325130D	•	•
15	•	•	EE325150D	•	•
<b>Piattaforma Platform</b>	•	•	• Ø 3,4	•	•

### NT

Connessione esagono esterno  
External hexagon connection



	Ø 3,25	Ø 4	Ø 5
H			
8,5	EE325085T	EE400085T	EE500085T
10	EE325100T	EE400100T	EE500100T
11,5	EE325115T	EE400115T	EE500115T
13	EE325130T	EE400130T	EE500130T
15	EE325150T	EE400150T	

Piattaforma  
Platform



Ø 3,4



Ø 4,1



Ø 5

### NT Compatibile Certain

Connessione esagono interno  
Internal hexagon connection



	Ø 3,25	Ø 4	Ø 5
H			
8,5	EI325085T	EI400085T	EI500085T
10	EI325100T	EI400100T	EI500100T
11,5	EI325115T	EI400115T	EI500115T
13	EI325130T	EI400130T	EI500130T
15	EI325150T	EI400150T	

Piattaforma  
Platform



Ø 3,4



Ø 4,1



Ø 5

## IMPIANTI TAPERED ACTION

COMPATIBLE NOBEL ACTIVE® Ø 3,5 NP (DAL Ø 3,5 AL Ø 5,2 COMPRESI)  
 COMPATIBLE NOBEL ACTIVE® DIAMETER 3,5 NP (FROM Ø 3,5 TO Ø 5,2)

DOPPIA connessione  
 CONICA esagonale  
 DOUBLE hexagonal  
 connection CONIC



Connessione  
 Connection

	Ø 3,1	Ø 3,5	Ø 3,8	Ø 4,2	Ø 4,7	Ø 5,2
H						
6	•	•	•	EI420060P	EI470060P	EI520060P
7,5	•	EI350080P	EI380080P	EI420075P	EI470075P	EI520075P
9	EI310090P	EI350090P	EI380090P	EI420090P	EI470090P	EI520090P
11	EI310110P	EI350110P	EI380110P	EI420110P	EI470110P	EI520110P
13	EI310130P	EI350130P	EI380130P	EI420130P	EI470130P	EI520130P
15	EI310150P	EI350150P	EI380150P	EI420150P	EI470150P	
17	EI310170P					
	■ Ø 2,9	■ Ø 3,5	■ Ø 3,8	■ Ø 4,2	■ Ø 4,7	■ Ø 5,2

## IMPIANTI AGO NEEDLE IMPLANTS



	Ø 2,4
H	
10	AG240100I
11,5	AG240115I
13	AG240130I
15	AG240150I

## IMPIANTI COMPATIBILI ASTRA-TECH™ COMPATIBLE ASTRA-TECH™ Ø 3,5



	Ø 3	Ø 3,4	Ø 3,8	Ø 4,5	Ø 5,5
H					
8	•	EI340080L	EI380080L	EI450080L	EI550080L
9,5	•	EI340095L	EI380095L	EI450095L	EI550095L
11	•	EI340110L	EI380110L	EI450110L	EI550110L
13	•	EI340130L	EI380130L	EI450130L	EI550130L
15	•	EI340150L	EI380150L	EI450150L	•
Connessione Connection	•	■ Ø 3,4	■ Ø 3,8	■ Ø 4,5	■ Ø 5,5







STRADA BATTAGLIA 127  
35020 ALBIGNASEGO (PD)  
ITALY

[info@rmedical.it](mailto:info@rmedical.it)  
[www.rmedical.it](http://www.rmedical.it)